

2020.09.15 - Oxygen Server Microscope 5.5.7 / Oxygen Immersion 5.5.6

Here you can find information about what is new in Oxygen Server Microscope 5.5.7 and Oxygen Immersion version 5.5.6.

On this page:


- 1 [Shared Features](#)
 - 1.1 [Camera FPS Values](#)
 - 1.2 [Fixed Problems and Improvements](#)

Shared Features

The functionality described in this section is available both in Oxygen Server Microscope and Oxygen Immersion.

Camera FPS Values

Now on the **Camera** panel, the camera FPS values are displayed for hardware (hw) and draw:



Fixed Problems and Improvements

The following fixes for the known problems and improvements have been implemented:

- Some fixes have been done resulting in system stability and speed.
- Previously the system had an expiration period starting from the distributive creation date, after this period the system could not be started. Now, this limitation is removed.