

Oxygen Immersion 3.16.15.1381 release notes

New release can be found here [Dropbox\Download\MboxIG\InstallOxygenImmersion3.16.15.1381Beta.exe](#)

What's new in this version:

- Rx endswitch monitor added. IMPORTANT: You need to update MicroscopeDrivers.ini. Details can be found [here](#).
- Parallel photo sets saving (gives ~30% scanning procedure speedup)
- Fixed Force mono checkbox: in previous versions even with Force mono ON left camera could significantly slow down scanning process
- Detailed warning message in case of very long exposure that can slow down scanning