

Alignment

Alignment is configuring of the MBox hardware functioning in a bundle with the corresponding software.

Some functions related to the MBox alignment can be accessed via the **Alignment** menu.



To use the **Alignment** menu, all the conditions should be met:

- You should use the system under the "Developer" HASP key or HASP key for Administrators (HASP key must have feature 108)
- The system should be launched by the shortcut of oxygen.exe with the "/alignmentmenu" flag

Name	Date modified	Type	Size
 Oxygen Immersion 5.1.9.3194 Beta	21.05.2019 12:55	Shortcut	2 KB
 Oxygen Server Microscope (with Alignment menu)	21.05.2019 16:28	Shortcut	2 KB
 Oxygen Server Microscope 5.1.11.3203	21.05.2019 16:28	Shortcut	2 KB
 Oxygen Web Site	21.05.2019 16:28	Internet Shortcut	1 KB

Oxygen Server Microscope (with Alignment menu) Pro...

Security

Details

Previous Versions

General

Shortcut

Compatibility

 Oxygen Server Microscope (with Alignment menu)

Target type: Application

Target location: Oxygen Server Microscope 5.1.11.3203

Target: scope 5.1.11.3203\oxygen.exe" /alignmentmenu

Child pages:

- Supported Graphic Cards
- Network Configuration and Settings
- Workstation Configuration — June 2014 (MBoxIG)